

APPLICANT

Institution : CRCM IPC Other :

Research team/ Company :

Project Manager : ☎ ✉

Applicant : ☎ ✉

Approval number DRC/DRT :

Project Name :

Financement :

SLIDES FOR SCANNING

Number of slides to scan : Magnification : x20 x40

Excel list of slides has to be provided under penalty of a service not performed

Slides to be saved on: J:\scanalam Y:\ICEP Team's hard drive (to be provided with the slides to be scanned)

Slides Format : .ndpi .wfml (Calopix)

Histo-Paca publication in Workshop :

Anonymization

INSTRUCTIONS TO FOLLOW

- In order to be scanned, your slides must meet the following criteria:
 - The label and coverslip must not extend beyond the glass slide
 - The sections must be centered on the slide; any section or part of a section too close to the edges of the slide cannot be scanned!
 - The slides must be cleaned (to avoid a out of focus image!).
- Sections or areas to be scanned should be **circled when the tissue is not clearly visible**.
- Your glass slides must be arranged in the order of your excel list.
- The images will be stored for 6 months on ICEP servers (J:\ or Y:\) and will then be deleted, so remember to make a copy!
- The glass slides must be retrieved within **15 days following the scanning of the slides**

ICEP

Date of scanning :

Save Path :

FLUOROCHROMES USED

Number of markers :

DAPI : Yes No

Fluorochromes :

Filters to use : DAPI FITC PE APC

To be completed by the applicant			Reserved for ICEP	
Filter	Expected labeling (cell localization)	Expected marking intensity (low, medium, high)	Exposure time (ms)	gain
DAPI				
FITC				
PE				
APC				

Scan profile name :

Comments :